


RNAscope Protocol ver.1

Pre-Assay Tissue Preparation

1. Slicing fresh frozen tissue
 - a. Freeze specimen with dry ice or liquid nitrogen
 - i. (adaptation: using isopentane)
 - b. Cut 10-20 μM sections and mount onto Superfrost Plus Slides
 - i. (adaptation: cutting at 16 μM , using StellarSci Slides)
 - c. Keep sections at -20°C to dry
 - d. Store the slide box (wrapped air tight) **in a ziplock bag at -80°C**
 - i. Sections can be stored for 3 months prior to use
(Optional stop point)
2. Pre-assay tissue processing
 - a. Reagents needed: 4% NBF in 1x PBS (chilled), 200 mL 50% Etoh (RT), 200 mL 70% Etoh (RT), 600mL 100% Etoh (RT) – Etoh made with 1X PBS
 - b. Pre-chill 200mL of 10% Normal buffered formalin or 4% PFA (in 1x PBS) in a tissue tek staining dish at 4°C
 - i. (adaptation: using 10% Neutral Buffered Formalin)
 1. Contains methanol which helps to permeabilize cells
 - c. Remove slides from -80°C , place in slide holder, and immediately **immerse in staining dish with 4% PFA in 1X PBS for 15min at 4°C .**
 - d. Place slide holder in staining dish with **50% Etoh at RT for 5min**
 - e. Place slide holder in staining dish with **70% Etoh at RT for 5min**
 - f. Place slide holder in staining dish with **100% Etoh at RT for 5min**
 - g. (Repeat 2f) Perform another incubation in **100% Etoh at RT for 5 min.**
 - h. Remove slides and place on absorbent paper (section facing up) **and air dry for 5 min.**
 - i. Draw a barrier with the ImmEdge hydrophobic barrier pen. (0.75" x 0.75") – use template on next page (1 template per slide)
 - i. **place image below underneath slide)**



- j. Let the barrier dry completely ~1 min.
3. Pre-treatment processing
 - a. **Turn on hybridization oven and set to 40°C** (also place a sheet of humidifying paper on the control tray and wet with distilled water).
 - b. At RT, place the dried slides on the HybEZ slide rack, add **2-4 drops of pretreat 4 per section and incubate for 30 min at RT.**
 - i. NOTE: this step may need to be adjusted between 15-30 min. depending on assay starting at 30 min
 - ii. for chromogenic assay, this step is preceded by pretreat 1 for 10 minutes followed by a 1x rinse in dH₂O

- c. Remove liquid (by tapping or flicking) and immediately place slides in a tissue tek slide rack in a tissue tek staining dish filled with 1x PBS (move up and down 3-5 times)
 - d. (Repeat 3c) – **PROCEED IMMEDIATELY TO RNAscope ASSAY**
 - i. Slides can stay in 1x PBS for upto 15 min.
- (--Do not apply Pretreat 1 or Pretreat 2 to fresh, frozen sections when using the Fluorescent Multiplex Kit. Other sample types may need these treatments.)

4. *Fluorescent RNAscope assay*

- a. Notes: make sure to run positive (Ppb) and negative control probe (DapB) on target tissue sample with each assay (1 slide each)
- b. Prep chemicals and probes
 - i. Wash Buffer
 - 1. 3L of wash buffer
 - a. Warm 50x wash buffer to 40°C for 10-20 min
 - b. Mix 1 bottle of 50x was buffer (60mL) to 2.94L distilled water
 - ii. Probes
 - 1. Warm probes (CRF C1, EGFP C3, Ppb, and DapB) **to 40°C for 10 min, then cool to RT**
 - 2. Spin down C2 or C3 probes
 - 3. Mix probes at 1:1:50 (C2 or C3 to C1) – (may need to cut top of C1 probe to ensure correct amount)
 - a. Mixed probes can be stored at 4°C for up to 6 mos.
 - iii. Bring AMP1-4FL to room temp
 - iv. Hybridize probe
 - 1. Tap to remove excess liquid and add **~4 drops of appropriate probe** to slides
 - 2. Place in HybEZ oven tray and secure in oven. **Incubate for 2 HRS at 40°C** (ensure turn lock on oven is secure to prevent evap)
 - 3. Remove slides 1 AT A TIME, tap to remove liquid, and immediately immerse in 1X wash buffer for 2 min at RT. Agitate slides by moving up and down.
 - 4. Repeat step 3
 - v. Hybridize AMP1-FL
 - 1. Tap to remove excess liquid and add **~4 drops of AMP1-FL** to slides
 - 2. Place in HybEZ oven tray and secure in oven. **Incubate for 30 min at 40°C** (ensure turn lock on oven is secure to prevent evap)
 - 3. Remove slides 1 AT A TIME, tap to remove liquid, and immediately immerse in **1X wash buffer for 2 min at RT.** **Agitate** slides by moving up and down.
 - 4. Repeat step 3
 - vi. Hybridize AMP2-FL
 - 1. Tap to remove excess liquid and add **~4 drops of AMP2-FL** to slides

2. Place in HybEZ oven tray and secure in oven. **Incubate for 15 min at 40°C** (ensure turn lock on oven is secure to prevent evap)
 3. Remove slides 1 AT A TIME, tap to remove liquid, and immediately immerse in **1X wash buffer for 2 min at RT**. **Agitate** slides by moving up and down.
 4. Repeat step 3
- vii. Hybridize AMP3-FL
1. Tap to remove excess liquid and add **~4 drops of AMP3-FL** to slides
 2. Place in HybEZ oven tray and secure in oven. **Incubate for 30 min at 40°C** (ensure turn lock on oven is secure to prevent evap)
 3. Remove slides 1 AT A TIME, tap to remove liquid, and immediately immerse in **1X wash buffer for 2 min at RT**. **Agitate** slides by moving up and down.
 4. Repeat step 3
- viii. Hybridize AMP4-FL
1. Tap to remove excess liquid and add **~4 drops of AMP4 ALT B-FL** to slides
 - a. Alt-B has C1 at Atto 550nm and C3 at Atto 647nm excitation
 2. Place in HybEZ oven tray and secure in oven. **Incubate for 15 min at 40°C** (ensure turn lock on oven is secure to prevent evap)
 3. Remove slides 1 AT A TIME, tap to remove liquid, and immediately immerse in **1X wash buffer for 2 min at RT**. **Agitate** slides by moving up and down.
 4. Repeat step 3
- ix. Counterstain and Mount
1. Tap to remove excess liquid and add **~4 drops of DAPI** to slides
 2. Place in HybEZ oven tray and secure in oven. **Incubate for 30 sec at RT**
 3. Remove DAPI from slides and immediately add Prolong Gold Antifade (without Dapi) and coverslip
 4. **Store slides in dark at 4°C**